



Quantitative determination of coating thickness on ABS parts

EDXRF Analyzer: **EX-6600**

Equipped with Silicon Drift Detector (SDD)



Abstract

Three plastic parts made of Acrylonitrile Butadiene Styrene (ABS) and coated by multi layers coating consisting of: Cu, Ni and Cr, were analysed on Xenometrix EX-6600 system, to determine the thickness of the different coating layers.

Determination of the thickness of the different coating layers was done using a Fundamental Parameter method that does not require a set of calibration standards.

Objective

To establish a method for measuring the thickness of multi layers coating on ABS, in the absence of calibration standards.

Background

EDXRF is a fast and non-destructive technique that can quantify elemental contents in any type of sample such as solid, powder, or liquids or determine thickness of elemental layers on different samples. A few advantages of Xenometrix EDXRF analyzers are: 1) no sample preparation required and 2) non-skilled operators can perform the analysis thanks to the user friendly graphical user interface. An additional advantage is the Fundamental parameter software for quantitative analysis without the need of costly and not always available calibration standards for the analysis.

Analytical Configuration

Table 1: Instrumental analytical configuration

Instrument	EX-6600
Anode	Rh-Anode X-ray Tube, 60kV,300W
Detector	Silicon Drift Detector (SDD)
Environment	Air
Excitation mode	Direct excitation
Type of analysis	Quantitative analysis
Analysis time	100 sec
Sample Preparation	No sample preparation

Experimental and Results

Three coated plastic parts were received for measuring the thickness of the different coated layers. The coating is consisting of the following layers:

1. 100% Cu – Inner layer; the closest to the ABS
2. 100% Ni – middle layer
3. 100% Cr – Outer layer


Figure 1: Sketch of Cu, Ni and Cr layers on the ABS samples

100% Cr
100% Ni
100% Cu
ABS (substrate)

Each plastic part was placed “as is” in the sample tray and analyzed. Typical spectrum of part 1 is shown in figure 2.

Static precision tests were performed on the three samples to show the repeatability of the method and the performance of EX-6600 Analyzer. The precision was performed by acquiring the spectra five times without moving the sample between acquisitions. The precision results; individual results and measured average concentration ± 1 standard deviation, are shown in table #2, 3 and 4.

Figure 2: Typical spectrum of Cu, Ni and Cr multi layers coating on ABS (part 1)



APPLICATION NOTE # XE-2014-3247**Quantitative determination of coating thickness on ABS parts****Table 2:** Sample 1 precision test results

	Cu [μm]	Ni [μm]	Cr [μm]
Sample 1_A	14.55	10.36	0.09
Sample 1_B	13.11	10.15	0.09
Sample 1_C	13.67	10.25	0.09
Sample 1_D	14.70	10.51	0.1
Sample 1_E	12.24	9.88	0.09
Average \pm standard deviation	13.66 \pm 1.03	10.23 \pm 0.24	0.09 \pm 0.00

Table 3: Sample 2 precision test results


	Cu [μm]	Ni [μm]	Cr [μm]
Sample 2_A	9.07	16.68	0.47
Sample 2_B	9.79	17.11	0.47
Sample 2_C	9.69	17.01	0.47
Sample 2_D	9.79	17.12	0.46
Sample 2_E	9.39	16.94	0.46
Average \pm standard deviation	9.54 \pm 0.31	16.97 \pm 0.18	0.47 \pm 0.00

Table 4: Sample 3 precision test results

	Cu [μm]	Ni [μm]	Cr [μm]
Sample 3_A	14.61	18.21	0.37
Sample 3_B	15.25	18.40	0.37
Sample 3_C	14.95	18.44	0.37
Sample 3_D	15.94	18.65	0.37
Sample 3_E	15.78	18.55	0.37
Average \pm standard deviation	15.31 \pm 0.56	18.45 \pm 0.17	0.37 \pm 0.00

Quantitative determination of coating thickness on ABS parts

Figure 3: Comparison between sample 1 (main red spectrum) and sample 2 (blue contour spectrum)



Spectra comparison shown in figure 3 above strongly suggests that the Cu layer in sample 1 is thicker than sample 2.

CONCLUSION

This application report shows how simple and rapid the determination of multi layer coating thickness can be performed on Xenometrix EDXRF analyzer, model EX-6600 SDD combined with Fundamental Parameter calculations. In the absence of calibration standards, fundamental parameter technique is the only solution for EDXRF quantitative analysis.